Centre Number

Candidate Number

Name

UNIVERSITY OF CAMBRIDGE INTERNATIONAL EXAMINATIONS General Certificate of Education Advanced Subsidiary Level and Advanced Level

CHEMISTRY 9701/03

Paper 3 Practical Test

May/June 2006

1 hour 15 minutes

Candidates answer on the Question Paper.

Additional Materials: As listed in the Instructions to Supervisors.

READ THESE INSTRUCTIONS FIRST

Write your details, including practical session and laboratory where appropriate, in the boxes provided. Write in dark blue or black pen.

You may use a soft pencil for any diagrams, graphs or rough working.

Do not use staples, paper clips, highlighters, glue or correction fluid.

Answer all questions.

You are advised to show all working in calculations.

Use of a Data Booklet is unnecessary.

Tests for gases are provided on page 8.

You may use a calculator.

At the end of the examination, fasten all your work securely together.

The number of marks is given in brackets [] at the end of each question or part question.

LABOR	LABORATORY	
For Exam	iner's Use	
I OI LAAIII	11161 3 036	
1		
•		
2		
TOTAL		

SESSION

This document consists of **7** printed pages and **1** blank page.

SP (SJF3983) T1246 © UCLES 2006



[Turn over

Examiner's Use

1 FA 1 is an aqueous bleach named 'Superclean'. FA 2 is a solution containing 19.78 g dm⁻³ of sodium thiosulphate, Na₂S₂O₃.

You are also provided with: aqueous potassium iodide, KI, $1.0\,\mathrm{mol}\,\mathrm{dm}^{-3}$ sulphuric acid, $\mathrm{H}_2\mathrm{SO}_4$, starch indicator solution, distilled water.

Many commercial bleaches contain compounds of chlorine as the 'active ingredient'. They are able to displace iodine from aqueous potassium iodide. This iodine can be titrated against sodium thiosulphate to determine the 'free chlorine' concentration for the bleach. You are to determine the 'free chlorine' concentration in 'Superclean'.

(a) Dilution of 'Superclean'

Pipette $25.0\,\mathrm{cm^3}$ of **FA 1** into the $250\,\mathrm{cm^3}$ graduated flask, labelled **FA 3**. Make the solution up to $250\,\mathrm{cm^3}$ with distilled water and mix the solution thoroughly. This solution is **FA 3**.

(b) Titration of iodine produced on reacting FA 3 with an excess of iodide ions

Rinse the 25 cm³ pipette used in **(a)**, firstly with water and then with **FA 3**. Pipette 25.0 cm³ of **FA 3** into a conical flask and, using a measuring cylinder, add to the flask 10 cm³ of aqueous potassium iodide and 10 cm³ of 1.0 mol dm⁻³ sulphuric acid.

Fill the burette with **FA 2** and titrate the displaced iodine until the colour of the solution becomes pale yellow. Then add 1 cm³ of starch indicator solution. Swirl the flask to mix the indicator and solution, then continue the titration until the blue colour of the starch/iodine complex disappears, leaving a colourless solution. This is the end-point.

Repeat the titration as many times as you think necessary to obtain accurate results.

Make certain that the recorded results show the precision of your practical work.

Table 1.1 Titration of liberated iodine with FA 2

final burette reading/cm ³		
initial burette reading/cm ³		
volume of FA 2 used/cm ³		

[7]

Summary

Show which results you used to obtain this volume of **FA 2** by placing a tick (\checkmark) under the readings in Table 1.1.



For Examiner's Use

You are advised to show full working in all parts of the calculations.

(c) Calculate how many moles of sodium thiosulphate, $Na_2S_2O_3$, were run from the burette during the titration. [A_r : Na, 23.0; S, 32.1; O, 16.0]

[2]

(d) Calculate how many moles of iodine, I₂, were present in the conical flask after **FA 3** reacted with an excess of iodide ions.

$$2S_2O_3^{2-}(aq) + I_2(aq) \longrightarrow S_4O_6^{2-}(aq) + 2I^{-}(aq)$$

[1]

(e) 1 mol of iodine is displaced by 1 mol of 'free chlorine'.

$$Cl_2(aq) + 2I^-(aq) \rightarrow 2Cl^-(aq) + I_2(aq)$$

Calculate the 'free chlorine' concentration, in mol dm⁻³, of **FA1**, the bleach 'Superclean'.

[3]



Examiner's
Use

(f) The *'free chlorine'* concentration of a second bleach, 'Germfree', is found by the same method.

 $50.0\,\mathrm{cm^3}$ of 'Germfree' is pipetted into a $250\,\mathrm{cm^3}$ graduated flask and the solution made up to $250\,\mathrm{cm^3}$ with distilled water. The 'free chlorine' concentration in the diluted solution is found to be $0.082\,\mathrm{mol\,dm^{-3}}$.

Calculate the 'free chlorine' concentration, in mol dm⁻³, in 'Germfree'.

[1]

(g) The cost of each bleach is as follows.

'Superclean' \$2.80 dm⁻³

'Germfree' \$1.80 dm⁻³

Use the data above and your answers to (e) and (f) to show which of the bleaches is better value for money.

[1]

[Total: 15]



BLANK PAGE



For Examiner's Use

2 You are provided with four aqueous solutions, FA 4, FA 5, FA 6 and FA 7.

Each solution contains one of the following.

an alcohol an aldehyde a carboxylic acid a ketone

You are to perform the tests below and from the results establish which type of organic compound is contained in each of **FA 4**, **FA 5**, **FA 6** and **FA 7**.

After each test discard the contents of the tubes into the 250 cm³ beaker, labelled organic waste. Rinse and re-use the tubes for the remaining tests.

Tests for gases are printed on page 8.

Record your results in the table on page 7. Where no reaction has taken place, write 'no change' in the appropriate box in the table.

test (a)	Place 1 cm depth of each of the solutions FA 4 , FA 5 , FA 6 and FA 7 into separate test-tubes. To each tube add a small quantity of magnesium powder or turnings. Identify any gas given off and record the test you used to make the identification.
test (b)	Place 1 cm depth of each of the solutions FA 4 , FA 5 , FA 6 and FA 7 into separate test-tubes. To each tube add a small quantity of powdered sodium carbonate. Identify any gas given off and record the test you used to make the identification.
test (c)	Place 1 cm depth of each of the solutions FA 4 , FA 5 , FA 6 and FA 7 into separate test-tubes. To each tube add 1 cm depth of 2,4-dinitrophenylhydrazine reagent.
test (d)	Place 1 cm depth of each of the solutions FA 4, FA 5, FA 6 and FA 7 into separate test-tubes. Place 2 cm depth of aqueous silver nitrate in a boiling-tube and add to it 1 cm depth of aqueous sodium hydroxide. This will produce a precipitate of silver oxide. Use a dropping pipette to add dilute aqueous ammonia to this mixture until the precipitate of silver oxide just dissolves. Do not add an excess of aqueous ammonia. To each of the tubes containing FA 4, FA 5, FA 6 and FA 7 add 1 cm depth of the silver-containing solution you have just prepared.
test (e)	Place 1 cm depth of each of the solutions FA 4 , FA 5 , FA 6 and FA 7 into separate boiling-tubes. To each tube add a few drops of acidified potassium dichromate(VI) to give a yellow-orange solution. Warm the tube gently.



For Examiner's Use

test	FA 4	FA 5	FA 6	FA 7
(a)				
(b)				
(c)				
(d)				
(e)				

[8]

Identify the type of organic compound present in each of the solutions FA 4, FA 5, FA 6 and FA 7 and complete the table below.

	type of organic compound contained in the solution	confirmed by the observations in test(s)
FA 4		
FA 5		
FA 6		
FA 7		

[2]

[Total: 10]



Tests for gases

gas	test and test result
ammonia, NH ₃	turns damp red litmus paper blue
carbon dioxide, CO ₂	gives a white ppt. with limewater (ppt. dissolves with excess CO ₂)
chlorine, Cl ₂	bleaches damp litmus paper
hydrogen, H ₂	'pops' with a lighted splint
oxygen, O ₂	relights a glowing splint
sulphur dioxide, SO ₂	turns potassium dichromate(VI) (aq) from orange to green

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

University of Cambridge International Examinations is part of the University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.

